

Apr 10, 2019

Working

Microscope video recording

Daniel Vaultot¹¹Station Biologique, Roscoff, France[dx.doi.org/10.17504/protocols.io.k24cygw](https://doi.org/10.17504/protocols.io.k24cygw)[Ecology of Marine Plankton \(ECOMAP\) team - Roscoff](#) [Roscoff Culture Collection](#)

Daniel Vaultot

Station Biologique, Roscoff, France



ABSTRACT

Protocol to record videos on SBR Microscope

GUIDELINES

Equipment needed

Hardware

- Inverted microscope IX-70 or IX-71
- Connecting tube between Microscope and C-mount
- [Infinity X Camera](#) (acquired in 2016)
- [Infinity 3-6 UR](#) (acquired in 2018)



Softwares

- [Image Capture software](#)
- [Video de Luxe software](#)

BEFORE STARTING

- Connect camera to microscope using C-mount connector
- Install Image Capture software on PC
- Connect USB cable to PC




First steps

- 1
 - Prepare your sample
 - Cultures flask
 - Drop between slide and coverslip and invert the slide so that the objective is on side of coverslip
 - Use x 40 objective for most cultures and x20 for dinoflagellates
 - Install and start [Infinity Capture software](#)


Acquisition with Infinity Capture

Acquisition with Infinity Capture


- 2 Advantage : Very easy to do
Disadvantage : Files are store as avi and are very very large...
- 3
 - Infinity X - Full resolution **1920x1080** (16:9) - Note : very large AVI files ! 30 sec = 700 Mo
 - Infinity 3UR-6 - Full resolution **2752 x 2192**




- 4
 - Infinity X - **800x608** (4:3) resolution for color videos - Note that only part of the field is visible but files are smaller than with full resolution
 - Infinity 3-6 UR **1024x768** ou **1600x1200**



- 5
 - Use 960x728 (4:3) resolution (binning) for BW videos - Note that the full field is visible
 -




- 6
 - Modify the AVI setups
 - Launch acquisition (usually only 10-20 sec at a time)



- 7
 - Need to buy the software (40 euros only)
 - Advantage : Much smaller files and the software allow to finalize the films

- 8
 - Setup the resolution in Image Capture
 - Shut down Image Capture
 - Launch Video de Luxe
 - Select in menu : Video Recording (Enregistrement video)
 - Slect the USB Camera - DShow (see next step)

- 9
 - Infinity X - Full resolution **1920x1080** (16:9)




- 10
 - **1280x720** (16:9) for selecting only center of camera field



11 Set up configuration as follows (item 3 in previous step)

- **Infinity X** - Resolution 1280 x720
- **Infinity 3-6UR** - Resolution 1024x768 (4:3)
- Ratio 16:9



Make the final movie with Video de Luxe


- 12
 - Create a new movie
 - Select the parts you want to add
 - Cut the scenes
 - Optimum duration is 20-30 sec max



- 13
 - Add titles and logos




- 14
 - Export to wmv format (Windows Movie - this is the most compact format without loss of resolution)
 - Choose a format with same aspect ratio than original
 - Do not add sound unless you want to add music
 - For 16:9 : 640 x 360 or 1280 x 720
 - For 4:3 : 640 x 480



Upload to YouTube

15



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited